


Anti-Secisbp2 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q96T21

Item Number: CM0027654



Introduction

Rabbit polyclonal antibody against human SECISBP2 (SBP2) validated for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat. Essential for selenoprotein synthesis research.

[Learn More](#)

Parameter	Value
Product Name	SECISBP2 Rabbit pAb
Remarks/Alias	SBP2; THMA1; SECISBP2
Species	Human
GeneID (Human)	79048
GeneID	79048
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 585-854 of human SECISBP2 (NP_076982.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96T21
Protein Weight	95kDa
Shipping	Ice bag