


Anti-Sel1L Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Ubv2

Item Number: CM0016338



Introduction

Rabbit polyclonal antibody against human SEL1L for Western blot, IF/ICC, and ELISA. Reacts with human, mouse, and rat samples. Suitable for studies of ER quality control and Notch signaling.

[Learn More](#)

Parameter	Value
Product Name	SEL1L Rabbit pAb
Remarks/Alias	Hrd3; SEL1L1; PRO1063; SEL1-LIKE; SEL1L
Species	Human
Gene ID (Human)	6400
Gene ID	6400
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-190 of human SEL1L (NP_005056.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UBV2
Protein Weight	89kDa
Shipping	Ice bag