


Anti-Sele Polyclonal Antibody For Wb And Elisa - P16581

Item Number: CM0018918



Introduction

High-quality rabbit polyclonal antibody against human SELE (E-selectin), validated for WB and ELISA, cross-reacts with mouse and rat. Ideal for studying leukocyte adhesion and endothelial biology.

[Learn More](#)

Parameter	Value
Product Name	SELE Rabbit pAb
Remarks/Alias	ELAM; ESEL; CD62E; ELAM1; LECAM2; selectin-e; SELE
Species	Human
Gene ID	6401
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-230 of human SELE (NP_000441.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P16581
Protein Weight	67kDa
Shipping	Ice bag