

Anti-Selenbp1 Rabbit Polyclonal Antibody For Wb, Elisa - Q13228

Item Number: CM0023294

Introduction

Rabbit polyclonal antibody targeting human SELENBP1 (Methanethiol oxidase), suitable for Western blot and ELISA. Cross-reacts with mouse. Useful for selenium biology, methanethiol metabolism, and cancer biomarker research.

[Learn More](#)

Parameter	Value
Product Name	SELENBP1 Rabbit pAb
Remarks/Alias	MTO; LPSB; SP56; hSBP; EHMT0; SBP56; HEL-S-134P; SELENBP1
Species	Human
Gene ID	8991
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 243-472 of human SELENBP1 (NP_003935.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q13228
Protein Weight	52kDa
Shipping	Ice bag