

Anti-Selenom Rabbit Polyclonal Antibody For Wb And Elisa - Q8Wwx9

Item Number: CM0019881

Introduction

SELENOM rabbit polyclonal antibody for WB and ELISA. Detects human and mouse SELENOM, a 16kDa selenoprotein involved in disulfide bond formation in cytoplasm, ER, and Golgi. Research use.

[Learn More](#)

Parameter	Value
Product Name	SELENOM Rabbit pAb
Remarks/Alias	SELM; SEPM; SELENOM
Species	Human
GeneID (Human)	140606
GeneID	140606
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 50-130 of human SELENOM (NP_536355.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q8WWX9
Protein Weight	16kDa
Shipping	Ice bag