

Anti-Selenop/Selenoprotein P Polyclonal Antibody For Wb, If/Icc, Elisa - P49908

Item Number: CM0019024

Introduction

Rabbit polyclonal antibody against human SELENOP (Selenoprotein P), validated for WB, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Suitable for investigating selenium transport and antioxidant defense. Targets the 43 kDa secreted glycoprotein.

[Learn More](#)

Parameter	Value
Product Name	SELENOP Rabbit pAb
Remarks/Alias	SeP; SELP; SEPP; SEPP1; SELENOP
Species	Human
Gene ID (Human)	6414
Gene ID	6414
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human,Mouse,Rat
SWISS	P49908
Protein Weight	43KDa
Shipping	Ice bag