

Anti-Sema3B Rabbit Polyclonal Antibody For Wb And Elisa - Q13214

Item Number: CM0030357

Introduction

Rabbit polyclonal antibody targeting human SEMA3B, validated for Western blot and ELISA. Detects human and mouse semaphorin-3B, a secreted axon guidance cue. Ideal for neural development studies.

[Learn More](#)

Parameter	Value
Product Name	SEMA3B Rabbit pAb
Remarks/Alias	SemaA; SEMA5; SEMAA; semaV; LUCA-1; SEMA3B
Species	Human
GeneID (human)	7869
GeneID	7869
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 569-748 of human SEMA3B (NP_001005914.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q13214
Protein Weight	83kDa
Shipping	Ice bag