

Anti-Sema3C Polyclonal Antibody For Wb, If/Icc, Elisa - Q99985

Item Number: CM0015556

Introduction

SEMA3C rabbit polyclonal antibody validated for WB, IF/ICC, and ELISA. Detects human and mouse semaphorin-3C, a key regulator of cardiovascular and neural development. Suitable for developmental research.

[Learn More](#)

Parameter	Value
Product Name	SEMA3C Rabbit pAb
Remarks/Alias	SemE; SEMAE; SEMA3C
Species	Human
GeneID	10512
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-330 of human SEMA3C (NP_006370.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q99985
Protein Weight	85kDa
Shipping	Ice bag