

Anti-Sema4A (Semaphorin-4A) Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9H3S1

Item Number: CM0021093

Introduction

SEMA4A Rabbit Polyclonal Antibody validated for WB, IF/ICC, ELISA across human, mouse, rat. Ideal for studying semaphorin signaling in neurobiology and immune regulation research and IVD assay development.

[Learn More](#)

Parameter	Value
Product Name	SEMA4A Rabbit pAb
Remarks/Alias	RP35; SEMB; SEMAB; CORD10; SEMA4A
Species	Human
GeneID (Human)	64218
GeneID	64218
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 320-580 of human SEMA4A (NP_071762.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H3S1
Protein Weight	84kDa
Shipping	Ice bag