

## Anti-Sema4C Polyclonal Antibody For Wb, Elisa - Q9C0C4

Item Number: CM0025340



### Introduction

Rabbit polyclonal antibody targeting human and mouse SEMA4C (Semaphorin-4C), validated for WB and ELISA. SEMA4C is a PLXNB2 receptor involved in axon guidance and cell signaling. Ideal for neuroscience and developmental biology research.

[Learn More](#)

| Parameter        | Value                                                                                                        |
|------------------|--------------------------------------------------------------------------------------------------------------|
| Product Name     | SEMA4C Rabbit pAb                                                                                            |
| Remarks/Alias    | SEMAF; SEMAI; SEMACL1; M-SEMA-F; SEMA4C                                                                      |
| Species          | Human                                                                                                        |
| Gene ID (Human)  | 54910                                                                                                        |
| Gene ID          | 54910                                                                                                        |
| Immunogen        | Recombinant protein corresponding to a sequence containing amino acids 749-833 of human SEMA4C (NP_060259.4) |
| Source           | Rabbit                                                                                                       |
| Category         | Polyclonal Antibodies                                                                                        |
| Application      | WB, ELISA                                                                                                    |
| Cross Reactivity | Human, Mouse                                                                                                 |
| SWISS            | Q9C0C4                                                                                                       |
| Protein Weight   | 93kDa                                                                                                        |
| Shipping         | Ice bag                                                                                                      |