

Anti-Sema4F Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - O95754

Item Number: CM0015065

Introduction

High-quality rabbit polyclonal antibody targeting human SEMA4F, suitable for Western blot, IHC-P, and ELISA. Recognizes mouse and rat orthologs. Ideal for neurogenesis and cell migration research.

[Learn More](#)

Parameter	Value
Product Name	SEMA4F Rabbit pAb
Remarks/Alias	S4F; SEMAM; SEMAW; M-SEMA; PRO2353; m-Sema-M; SEMA4F
Species	Human
Gene ID	10505 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 280-420 of human SEMA4F (NP_004254.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95754
Protein Weight	84kDa
Shipping	Ice bag