


Anti-Sema7A Polyclonal Antibody For Wb, Elisa - 075326

Item Number: CM0015960



Introduction

Rabbit polyclonal antibody against human SEMA7A (Semaphorin-7A, CD108), validated for WB and ELISA. SEMA7A regulates integrin signaling, immune responses, and axon guidance. Ideal for inflammation and neurobiology research.

[Learn More](#)

Parameter	Value
Product Name	SEMA7A Rabbit pAb
Remarks/Alias	JMH, CD108, SEMAL, CDw108, PFIC11, SEMAK1, H-Sema-L, H-SEMA-K1, SEMA7A
Species	Human
Gene ID (Human)	8482
Gene ID	8482
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 45-280 of human SEMA7A (NP_003603.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	075326
Protein Weight	75kDa
Shipping	Ice bag