

Anti-Semaphorin 3A Rabbit Polyclonal Antibody For Wb, Elisa - Q14563

Item Number: CM0016808

Introduction

Rabbit polyclonal antibody targeting human Semaphorin 3A for western blotting and ELISA. Cross-reacts with mouse and rat. Useful for studies of neuronal growth cone guidance and olfactory development.

[Learn More](#)

Parameter	Value
Product Name	Semaphorin 3A Rabbit pAb
Remarks/Alias	HH16; SemD; COLL1; SEMA1; SEMAD; SEMAL; coll-1; Hsema-I; SEMAIII; Hsema-III; Semaphorin 3A
Species	Human
Gene ID (Human)	10371
Gene ID	10371
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 300-400 of human Semaphorin 3A (NP_006071.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q14563
Protein Weight	89kDa
Shipping	Ice bag