


Anti-Semaphorin 4D / Cd100 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q92854

Item Number: CM0003295



Introduction

High-affinity rabbit monoclonal antibody targeting human Semaphorin 4D/CD100. Validated for WB, IHC-P, ELISA applications. Ideal for studying immune regulation, neuronal development, and semaphorin-plexin signaling in diagnostic research.

[Learn More](#)

Parameter	Value
Product Name	Semaphorin 4D / CD100 Rabbit mAb
Remarks/Alias	A8; GR3; BB18; CD100; COLL4; SEMAJ; coll-4; C9orf164; M-sema-G
Species	Human
GeneID (Human)	10507
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 100-200 of human Semaphorin 4D / CD100 (NP_006369.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	Q92854
Protein Weight	96kDa
Shipping	Ice bag