


Anti-Semaphorin-4D/Cd100 Rabbit Monoclonal Antibody For IF/ICC, Fc, Elisa - Q92854

Item Number: CM0007583



Introduction

Recombinant rabbit monoclonal antibody against human Semaphorin-4D/CD100. Validated for IF/ICC, FC, ELISA. Targets an 82/96 kDa cell surface receptor involved in neurogenesis and immune response. Ideal for neuroscience and immunology research.

[Learn More](#)

Parameter	Value
Product Name	Semaphorin-4D/CD100 Rabbit mAb
Remarks/Alias	SEMA4D; C9orf164; CD100; M-sema-G; SEMA; coll-4; semaphorin-4D; Semaphorin-4D/CD100
Species	Human
GeneID (human)	10507
GeneID	10507
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-734 of human Semaphorin-4D/CD100 (NP_006369.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF/ICC, FC, ELISA
Cross Reactivity	Human
SWISS	Q92854
Protein Weight	82kDa/96kDa
Shipping	Ice bag