

Anti-Semg1 Polyclonal Antibody For Wb, Elisa - P04279

Item Number: CM0020786



Introduction

Rabbit polyclonal antibody against human Semenogelin-1 (SEMG1), suitable for Western blot and ELISA. Reacts with human species; immunogen spans amino acids 312-462. Essential for reproductive biology and fertility research.

[Learn More](#)

Parameter	Value
Product Name	SEMG1 Rabbit pAb
Remarks/Alias	SEMG1; CT103; SEMG; SGI; dj172H20.2; semenogelin-1
Species	Human
GeneID	6406
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 312-462 of human SEMG1 (NP_002998.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P04279
Protein Weight	45kDa/52kDa
Shipping	Ice bag