

Anti-Semg2 Rabbit Polyclonal Antibody For Wb, Elisa - Q02383

Item Number: CM0016381

Introduction

Rabbit polyclonal antibody against human Semenogelin-2 (SEMG2) for WB and ELISA. Detects human and mouse SEMG2. Ideal for reproductive biology and semen coagulation studies.

[Learn More](#)

Parameter	Description
Product Name	SEMG2 Rabbit pAb
Remarks/Alias	SGII; SEMG2
Species	Human
Gene ID (Human)	6407
Gene ID	6407
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 383-582 of human SEMG2 (NP_002999.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q02383
Protein Weight	65kDa
Shipping	Ice bag