


Anti-Senp1 Rabbit Monoclonal Antibody For Wb, Elisa - Q9P0U3

Item Number: CM0007943



Introduction

Anti-SEN1 rabbit monoclonal antibody validated for WB and ELISA. Targets human sentrin-specific protease 1 (SuPr-2, SENP1), a 73 kDa SUMO pathway hydrolase critical for protein deSUMOylation. Suitable for studying transcriptional regulation and RNA methylation.

[Learn More](#)

Parameter	Value
Product Name	SEN1 Rabbit mAb
Remarks/Alias	SuPr-2; SEN1
Species	Human
Gene ID (Human)	29843
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 400-501 of human SEN1 (Q9P0U3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9P0U3
Protein Weight	73kDa
Shipping	Ice bag