


Anti-Senp1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9P0U3

Item Number: CM0014745



Introduction

Rabbit polyclonal antibody against human sentrin-specific protease 1 (SEN1). Validated for WB, IHC-P, ELISA; cross-reacts with mouse and rat. Ideal for SUMOylation studies and deSUMOylase research.

[Learn More](#)

Parameter	Value
Product Name	SEN1 Rabbit pAb
Remarks/Alias	SuPr-2; SEN1
Species	Human
GeneID	29843 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 220-420 of human SEN1 (NP_001254523.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9P0U3
Protein Weight	73kDa
Shipping	Ice bag