


Anti-Senp2 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Hc62

Item Number: CM0023469



Introduction

Validated rabbit polyclonal antibody against human SENP2 for Western blot, IHC paraffin, and ELISA. Reacts with human, mouse, rat. Suitable for SUMO pathway, Wnt signaling, and adipogenesis research. For research use only.

[Learn More](#)

Parameter	Value
Product Name	SENP2 Rabbit pAb
Remarks/Alias	AXAM2, SMT3IP2, SENP2
Species	Human
Gene ID (Human)	59343
Gene ID	59343
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 366-589 of human SENP2 (NP_067640.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9HC62
Protein Weight	68kDa
Shipping	Ice bag