

Anti-Sepn1 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Nzv5

Item Number: CM0013548

Introduction

Rabbit polyclonal antibody targeting SEPN1 (Selenoprotein N) for Western blot and ELISA. Cross-reacts with human and mouse. Ideal for oxidative stress, calcium homeostasis, and muscle biology studies.

[Learn More](#)

Parameter	Value
Product Name	SEPN1 Rabbit pAb
Remarks/Alias	RSS; CFTD; SELN; CMYP3; MDRS1; RSMD1; SEPN1
Species	Human
Gene ID	57190
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 341-590 of human SEPN1 (NP_065184.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9NZV5
Protein Weight	66kDa
Shipping	Ice bag