

Anti-Sepsecs Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Hd40

Item Number: CM0028202



Introduction

SEPSECS rabbit polyclonal antibody for WB, IF/ICC, ELISA; detects human, mouse, rat SEPSECS, a selenocysteine synthase essential for selenoprotein biosynthesis and a known autoantigen in liver disease.

[Learn More](#)

Parameter	Value
Product Name	SEPSECS Rabbit pAb
Remarks/Alias	LP; SLA; SecS; PCH2D; SLA/LP; SLA-p35; SEPSECS
Species	Human
GeneID (Human)	51091
GeneID	51091
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 222-501 of human SEPSECS (NP_058651.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9HD40
Protein Weight	56kDa
Shipping	Ice bag