

Anti-Septin 11 Polyclonal Antibody For Wb, Elisa - Q9Nva2

Item Number: CM0018488



Introduction

Rabbit polyclonal antibody targeting human Septin-11, validated for WB and ELISA. Cross-reacts with mouse and rat. For cytoskeleton, cell division, and neuroscience studies.

[Learn More](#)

Parameter	Value
Product Name	Septin 11 Rabbit pAb
Remarks/Alias	SEPT11; Septin 11
Species	Human
GeneID (Human)	55752
GeneID	55752
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-120 of human Septin 11 (NP_060713.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NVA2
Protein Weight	49kDa
Shipping	Ice bag