

Anti-Septin 2 Polyclonal Antibody For Wb, If/Icc, Elisa - Q15019

Item Number: CM0026146



Introduction

High-quality anti-Septin 2 rabbit polyclonal antibody for Western blot, immunofluorescence, and ELISA. Reacts with human and mouse Septin-2. Recombinant protein immunogen. Ideal for cell division and cytoskeleton research.

[Learn More](#)

Parameter	Value
Product Name	Septin 2 Rabbit pAb
Remarks/Alias	DIFF6; NEDD5; SEPT2; NEDD-5; Pnutl3; hNedd5; Septin 2
Species	Human
GeneID (Human)	4735
GeneID	4735
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 112-361 of human Septin 2 (NP_004395.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q15019
Protein Weight	41kDa
Shipping	Ice bag