

Anti-Septin 2 Rabbit Monoclonal Antibody For Wb, Elisa - Q15019

Item Number: CM0011403



Introduction

Rabbit monoclonal antibody against human Septin 2 (NEDD-5), validated for WB and ELISA. Reacts with human, mouse, and rat. Septin-2 is a filament-forming GTPase involved in mitosis, ciliogenesis, and actin cytoskeleton organization. Ideal for cell biology research.

[Learn More](#)

Parameter	Specification
Product Name	Septin 2 Rabbit mAb
Remarks/Alias	DIFF6; NEDD5; SEPT2; NEDD-5; Pnutl3; hNedd5; Septin 2
Species	Human
Gene ID	4735
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 202-298 of human Septin 2 (Q15019)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15019
Protein Weight	41kDa
Shipping	Ice bag