

Anti-Septin 3 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Uh03

Item Number: CM0021226



Introduction

Septin 3 Rabbit Polyclonal Antibody (41kDa) validated for WB, IF/ICC, ELISA in Human, Mouse, Rat. Targets neuronal-specific septin-3, a brain-specific cytoskeletal GTPase involved in cytokinesis. Ideal for neurobiology research.

[Learn More](#)

Parameter	Value
Product Name	Septin 3 Rabbit pAb
Remarks/Alias	SEP3; SEPT3; bK250D10.3; Septin 3
Species	Human
GeneID	55964
Immunogen	Recombinant Protein: Recombinant fusion protein containing a sequence corresponding to amino acids 179-358 of human Septin 3 (NP_663786.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UH03
Protein Weight	41kDa
Shipping	Ice bag