


Anti-Septin 4 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - O43236

Item Number: CM0016864



Introduction

Rabbit polyclonal antibody against Septin 4 (SWISS: O43236), suitable for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with human, mouse, rat. Detects cytoskeletal GTPase involved in apoptosis and spermatogenesis.

[Learn More](#)

Parameter	Value
Product Name	Septin 4 Rabbit pAb
Remarks/Alias	H5, ARTS, MART, SEP4, CE5B3, SEPT4, PNUTL2, hucep-7, BRADEION, C17orf47, hCDCREL-2, Septin 4
Species	Human
Gene ID	5414
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43236
Protein Weight	55kDa
Shipping	Ice bag