

Anti-Septin 5 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q99719

Item Number: CM0015574



Introduction

Rabbit polyclonal Septin 5 antibody validated for WB, IHC-P, IF/ICC, and ELISA. Detects human, mouse, and rat Septin 5, a filament-forming cytoskeletal GTPase involved in cytokinesis and platelet secretion.

[Learn More](#)

Parameter	Value
Product Name	Septin 5 Rabbit pAb
Remarks/Alias	H5; SEPT5; CDCREL; PNUTL1; CDCREL1; CDCREL-1; HCDCREL-1; Septin 5
Species	Human
GeneID (Human)	5413
GeneID	5413
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-369 of human Septin 5 (NP_002679.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99719
Protein Weight	43kDa
Shipping	Ice bag