

Anti-Septin-6 Rabbit Polyclonal Antibody For Wb And Elisa - Q14141

Item Number: CM0027307



Introduction

Rabbit polyclonal antibody against human Septin-6 (Q14141), validated for WB and ELISA. Cross-reacts with mouse and rat. Ideal for cytoskeleton, cell division, and spermatogenesis studies.

[Learn More](#)

Parameter	Value
Product Name	Septin 6 Rabbit pAb
Remarks/Alias	SEP2; SEPT2; SEPT6; Septin 6
Species	Human
Gene ID	23157
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 350-434 of human Septin 6 (NP_055944.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14141
Protein Weight	50kDa
Shipping	Ice bag