

Anti-Septin 7 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q16181

Item Number: CM0026276



Introduction

Rabbit polyclonal antibody against human Septin 7 (Q16181), validated for WB, IHC-P, ELISA. Reacts with human, mouse, rat. Recombinant immunogen targeting amino acids 307-436. Suitable for cytoskeleton, mitosis, and ciliogenesis research.

[Learn More](#)

Parameter	Value
Product Name	Septin 7 Rabbit pAb
Remarks/Alias	CDC3, CDC10, SEPT7, SEPT7A, NBLA02942, Septin 7
Species	Human
GeneID (Human)	989
GeneID	989
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 307-436 of human Septin 7 (NP_001011553.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16181
Protein Weight	51kDa
Shipping	Ice bag