

Anti-Septin 8 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q92599

Item Number: CM0018283



Introduction

Rabbit polyclonal antibody targeting human Septin 8, suitable for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Ideal for studying cytoskeleton dynamics, vesicle trafficking, and Alzheimer's disease-related pathways.

[Learn More](#)

Parameter	Value
Product Name	Septin 8 Rabbit pAb
Remarks/Alias	SEP2; SEPT8; Septin 8
Species	Human
GeneID (human)	23176
GeneID	23176
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 348-427 of human Septin 8 (NP_001287728.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92599
Protein Weight	56kDa
Shipping	Ice bag