

Anti-Serca1/Atp2A1 Monoclonal Antibody For Wb, If-P, Elisa - O14983

Item Number: CM0003277

Introduction

Recombinant rabbit monoclonal antibody targeting human SERCA1/ATP2A1 (O14983). Suitable for Western blot, immunofluorescence, and ELISA. Cross-reacts with mouse and rat. Ideal for calcium pump research.

[Learn More](#)

Parameter	Value
Product Name	SERCA1/ATP2A1 Rabbit mAb
Remarks/Alias	ATP2A; SERCA1; SERCA1/ATP2A1
Species	Human
GeneID (Human)	487
GeneID	487
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 900-1001 of human SERCA1/ATP2A1 (O14983).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	O14983
Protein Weight	95kDa/110kDa
Shipping	Ice bag