

Anti-Serca1/Atp2A1 Polyclonal Antibody For Wb - O14983

Item Number: CM0012915

Introduction

Rabbit polyclonal antibody targeting SERCA1/ATP2A1. Validated for Western blot and ELISA, cross-reacts with mouse. Suitable for research on calcium transport in skeletal muscle.

[Learn More](#)

Parameter	Value
Product Name	SERCA1/ATP2A1 Rabbit pAb
Remarks/Alias	ATP2A; SERCA1; SERCA1/ATP2A1
Species	Human
GeneID (Human)	487
GeneID	487
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 901-1001 of human SERCA1/ATP2A1 (NP_775293.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	O14983
Protein Weight	110kDa
Shipping	Ice bag