

Anti-Serca2/Atp2A2 Monoclonal Antibody For Wb, If/Icc, Ip, Elisa - P16615

Item Number: CM0002974



Introduction

Rabbit monoclonal antibody against SERCA2/ATP2A2, validated for WB, IF/ICC, IP, ELISA; cross-reacts with human, mouse, rat. Ideal for studies of calcium transport and ER/SR function.

[Learn More](#)

Parameter	Value
Product Name	SERCA2/ATP2A2 Rabbit mAb
Remarks/Alias	DD, DAR, ATP2B, SERCA2
Species	Human
GeneID (Human)	488
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 943-1042 of human SERCA2/ATP2A2 (P16615).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P16615
Protein Weight	115 kDa
Shipping	Ice bag