

Anti-Serca2/Atp2A2 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P16615

Item Number: CM0014713



Introduction

Rabbit polyclonal antibody against human SERCA2/ATP2A2, validated for Western blot, IHC-P, IF/ICC, and ELISA; cross-reacts with mouse and rat. Suitable for calcium signaling and ER stress research.

[Learn More](#)

Parameter	Value
Product Name	SERCA2/ATP2A2 Rabbit pAb
Remarks/Alias	DD; DAR; ATP2B; SERCA2; SERCA2/ATP2A2
Species	Human
Gene ID (Human)	488
Gene ID	488
Immunogen	Recombinant fusion protein containing a sequence within amino acids 1-80 of human SERCA2/ATP2A2 (NP_733765.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P16615
Protein Weight	115kDa
Shipping	Ice bag