

Anti-Serca2/Atp2A2 Rabbit Polyclonal Antibody For Wb, Ihc-P - P16615

Item Number: CM0015811



Introduction

Rabbit polyclonal antibody targeting human SERCA2/ATP2A2 (P16615), validated for Western blot, immunohistochemistry (paraffin), and ELISA. Cross-reacts with mouse and rat. Ideal for calcium transport and ER/SR membrane studies.

[Learn More](#)

Parameter	Value
Product Name	SERCA2/ATP2A2 Rabbit pAb
Remarks/Alias	DD; DAR; ATP2B; SERCA2; SERCA2/ATP2A2
Species	Human
GeneID (human)	488
GeneID	488
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 111-253 of human SERCA2/ATP2A2 (NP_733765.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P16615
Protein Weight	115kDa
Shipping	Ice bag