

Anti-Serine Racemase Rabbit Monoclonal Antibody For Wb And Elisa - Q9Gzt4

Item Number: CM0003854



Introduction

Anti-Serine racemase rabbit monoclonal antibody, reactive with human, for WB and ELISA applications. Serine racemase converts L-serine to D-serine, a co-agonist of NMDA receptors; relevant in neurology and schizophrenia studies.

[Learn More](#)

Parameter	Value
Product Name	Serine racemase Rabbit mAb
Alias	ILV1, ISO1, Serine racemase
Species	Human
GeneID (Human)	63826
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of human Serine racemase (NP_068766.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9GZT4
Protein Weight	37kDa
Shipping	Ice bag