

Anti-Serine Racemase (Srr) Rabbit Polyclonal Antibody For Wb, Elisa - Q9Gzt4

Item Number: CM0025463



Introduction

Rabbit polyclonal antibody against human serine racemase (SRR) for WB and ELISA. Cross-reacts with human, mouse, rat. Molecular weight: 37kDa. Ideal for D-serine synthesis and NMDA receptor studies.

[Learn More](#)

Parameter	Value
Product Name	SRR Rabbit pAb
Remarks/Alias	ILV1; ISO1; SRR
Species	Human
GeneID (Human)	63826
GeneID	63826
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-340 of human SRR (NP_068766.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9GZT4
Protein Weight	37kDa
Shipping	Ice bag