

Anti-Serotonin Transporter Polyclonal Antibody For Wb, Ihc-P, Elisa - Q60857

Item Number: CM0030214

Introduction

This rabbit polyclonal antibody targets the N-terminal region of mouse serotonin transporter (SERT/Slc6a4). Validated in WB, IHC-P, ELISA; cross-reacts with mouse and rat. Ideal for serotonin signaling studies.

[Learn More](#)

Parameter	Value
Product Name	Serotonin transporter Rabbit pAb
Remarks/Alias	Htt; Sert; 5-HTT
Species	Mouse
GenelD (Human)	6532
GenelD	15567
Immunogen	Recombinant protein (Recombinant fusion protein containing a sequence corresponding to amino acids 1-85 of mouse Serotonin transporter (NP_034614.2))
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q60857
Protein Weight	70 kDa
Shipping	Ice bag