

Anti-Serpina10 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9UK55

Item Number: CM0026829



Introduction

High-quality anti-SERPINA10 rabbit polyclonal antibody validated for Western blot, IHC-P, and ELISA. Recognizes human, mouse, and rat SERPINA10 (Protein Z-dependent protease inhibitor), key for coagulation research.

[Learn More](#)

Parameter	Value
Product Name	SERPINA10 Rabbit pAb
Remarks/Alias	PZI; ZPI; SERPINA10
Species	Human
GenelD (human)	51156
GenelD	51156
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UK55
Protein Weight	51kDa
Shipping	Ice bag