

## Anti-Serpina12 Polyclonal Antibody For Wb, Elisa - Q8lw75

Item Number: CM0030720



### Introduction

Rabbit polyclonal antibody targeting human SERPINA12 (Vaspin). Validated for WB and ELISA; cross-reacts with human and mouse. Suitable for insulin action and adipose tissue studies.

[Learn More](#)

| Parameter        | Value                                                                                                                                     |
|------------------|-------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | SERPINA12 Rabbit pAb                                                                                                                      |
| Remarks/Alias    | OL-64; SERPINA12                                                                                                                          |
| Species          | Human                                                                                                                                     |
| GeneID           | 145264                                                                                                                                    |
| Immunogen        | Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 20-320 of human SERPINA12 (NP_776249.1) |
| Source           | Rabbit                                                                                                                                    |
| Category         | Polyclonal Antibodies                                                                                                                     |
| Application      | WB, ELISA                                                                                                                                 |
| Cross Reactivity | Human, Mouse                                                                                                                              |
| SWISS            | Q8IW75                                                                                                                                    |
| Protein Weight   | 47kDa                                                                                                                                     |
| Shipping         | Ice bag                                                                                                                                   |