

Anti-Serpina6 Rabbit Polyclonal Antibody For Wb, Elisa - P08185

Item Number: CM0028376



Introduction

High-quality rabbit polyclonal antibody against human SERPINA6 (corticosteroid-binding globulin/transcortin), validated for Western blot and ELISA. Ideal for investigating glucocorticoid transport and steroid hormone biology.

[Learn More](#)

Parameter	Value
Product Name	SERPINA6 Rabbit pAb
Remarks/Alias	CBG; SERPINA6
Species	Human
Gene ID (Human)	866
Gene ID	866
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 136-405 of human SERPINA6 (NP_001747.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P08185
Protein Weight	45kDa
Shipping	Ice bag