

Anti-Serpina9 Polyclonal Antibody For Wb, Elisa - Q86Wd7

Item Number: CM0020747

Introduction

High-quality rabbit polyclonal antibody against human SERPINA9 (Serpine A9), validated for Western blot and ELISA. Suitable for studying germinal center B-cell biology and serpin function.

[Learn More](#)

Parameter	Value
Product Name	SERPINA9 Rabbit pAb
Remarks/Alias	GCET1; SERPINA11; SERPINA11b; SERPINA9
Species	Human
GeneID (human)	327657
GeneID	327657
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 110-210 of human SERPINA9 (NP_783866.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q86WD7
Protein Weight	47kDa
Shipping	Ice bag