

Anti-Serpinb1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P30740

Item Number: CM0029909



Introduction

Anti-SERPINB1 rabbit polyclonal antibody for detection of leukocyte elastase inhibitor in WB, IHC-P, IF/ICC, ELISA. Cross-reacts with human, mouse, rat. Supports studies of innate immunity and protease inhibition.

[Learn More](#)

Parameter	Value
Product Name	SERPINB1 Rabbit pAb
Remarks/Alias	EI; LEI; PI2; MNEI; PI-2; HEL57; M/NEI; ELANH2; HEL-S-27; SERPINB1
Species	Human
GeneID (Human)	1992
GeneID	1992
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 100-379 of human SERPINB1 (NP_109591.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P30740
Protein Weight	43kDa
Shipping	Ice bag