

Anti-Serpinb13 Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Uiv8

Item Number: CM0017244

Introduction

Rabbit polyclonal antibody targeting SERPINB13 (Hurpin/Headpin). Validated for WB, IF/ICC, ELISA; reacts with human, mouse, rat. Predicted MW 44 kDa. Ideal for skin biology and keratinocyte differentiation studies.

[Learn More](#)

Parameter	Value
Product Name	SERPINB13 Rabbit pAb
Remarks/Alias	HUR7; PI13; headpin; HSHUR7SEQ; SERPINB13
Species	Human
GeneID (Human)	5275
GeneID	5275
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-140 of human SERPINB13 (NP_036529.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UIV8
Protein Weight	44kDa
Shipping	Ice bag