

Anti-Serpib2/Pai-2 Rabbit Polyclonal Antibody For Wb And Elisa - P05120

Item Number: CM0020520



Introduction

Rabbit polyclonal antibody against human SERPINB2 (PAI-2) for Western blot (WB) and ELISA. This reagent supports research into urokinase-type plasminogen activator inhibition, serpin function, and related IVD applications.

[Learn More](#)

Parameter	Value
Product Name	SERPINB2 Rabbit pAb
Remarks/Alias	PAI; PAI2; PAI-2; PLANH2; HsT1201; SERPINB2
Species	Human
Gene ID (Human)	5055
Gene ID	5055
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-415 of human SERPINB2 (NP_002566.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P05120
Protein Weight	47kDa
Shipping	Ice bag