

Anti-Serpinb6 Polyclonal Antibody For Wb, Elisa - P35237

Item Number: CM0029239

Introduction

SERPINB6 rabbit polyclonal antibody validated for WB and ELISA, reacting with human, mouse, and rat. Ideal for studying serine protease inhibition, hearing loss, and cytoplasmic protein regulation.

[Learn More](#)

Parameter	Value
Product Name	SERPINB6 Rabbit pAb
Remarks/Alias	CAP; PI6; PTI; PI-6; SPI3; DFNB91; MSTP057; SERPINB6
Species	Human
GeneID (Human)	5269
GeneID	5269
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 180-300 of human SERPINB6 (NP_004559.4).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P35237
Protein Weight	43kDa
Shipping	Ice bag