

Anti-Serpinb9 (Pi9/Cap-3) Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P50453

Item Number: CM0030048



Introduction

IVD-grade SERPINB9 rabbit polyclonal antibody validated for WB, IF/ICC, ELISA. Recognizes human, mouse, rat serpin B9 (PI9/CAP-3), a cytoplasmic granzyme B inhibitor. Ideal for apoptosis and immune regulation research.

[Learn More](#)

Parameter	Value
Product Name	SERPINB9 Rabbit pAb
Remarks/Alias	PI9; CAP3; PI-9; CAP-3; SERPINB9
Species	Human
Gene ID (Human)	5272
Gene ID	5272
Immunogen	Recombinant fusion protein containing amino acids 170-308 of human SERPINB9 (NP_004146.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P50453
Protein Weight	42kDa
Shipping	Ice bag