

Anti-Serpinc1 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P01008

Item Number: CM0003863



Introduction

High-quality rabbit monoclonal antibody specific for human SERPINC1 (Antithrombin-III). Validated for Western Blot, IHC-P, and ELISA; cross-reacts with human, mouse, and rat. Ideal for coagulation research.

[Learn More](#)

Parameter	Value
Product Name	SERPINC1 Rabbit mAb
Remarks/Alias	AT3, AT3D, ATIII, THPH7, ATIII-R2, ATIII-T1, ATIII-T2, SERPINC1
Species	Human
Gene ID (Human)	462
Gene ID	462
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human SERPINC1 (P01008).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P01008
Protein Weight	58kDa
Shipping	Ice bag