


Anti-Serpine2 Rabbit Polyclonal Antibody For Wb, Elisa - P07093

Item Number: CM0027906



Introduction

Rabbit polyclonal antibody recognizing human SERPINE2 (Glia-derived nexin, PN-1). Validated for WB and ELISA; cross-reacts with mouse. Ideal for studies of serine protease inhibition and neurogenesis.

[Learn More](#)

Parameter	Value
Product Name	SERPINE2 Rabbit pAb
Remarks/Alias	GDN; PI7; PN1; PNI; PI-7; PN-1; GDNPF; SERPINE2
Species	Human
GeneID (human)	5270
GeneID	5270
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-398 of human SERPINE2 (NP_006207.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	P07093
Protein Weight	44kDa
Shipping	Ice bag